

Aldasan® 2000 / Aldosan® 2000 - CH

Concentrate **Disinfecting and cleaning concentrate for surfaces, equipment, medical instruments and endoscopes**

For professional use

COMPOSITION

100 g solution contains:

9.8 g Formaldehyde, 9.9 g Glutaral, Sodium alkyl ether sulphate, Isotridecanol (ethoxylated)

RANGE OF APPLICATION

- **One-step disinfection and cleaning** of surfaces, equipment and instruments (medical devices included) in case of bacterial (incl. Tb), viral and fungal infections; also suitable for the disinfection of textiles like mops, cloths etc.
- Control of nosocomial infections.
- Disinfection and cleaning of floors and objects (e.g. bathroom mats and grids, shoes, socks etc.) contaminated by fungus.
- Disinfection and hygienic cleaning of all sanitary equipment.
- Can be used for ultra-sound disinfection.

MODE OF APPLICATION / DOSAGE

Aldasan® 2000/Aldosan2000-CH is a liquid concentrate with which any use-dilution required can be made by simply adding water (at max. room temperature). Always fill container with tap water first and then add disinfectant.

Efficacy	Concentration	Exposure time
Disinfection by wiping accord. to the VAH-procedures dated 2002 with mechanical effect, strong contamination		
bactericidal, levurocidal, fungicidal	1.00%	15 min
	0.75%	30 min
	0.50%	1 h
	0.25%	4 h
Instrument disinfection accord. to the VAH-procedures dated 2002, strong contamination		
bactericidal, levurocidal, fungicidal	2.00%	30 min
	1.00%	1 h
mycobactericidal	3.00%	30 min
	2.00%	1 h
Tests according to the DVV/RKI directive		
Virucidal against enveloped viruses (incl. HIV, HBV, HCV) According to the recommendations of RKI 1/2004*	0.50%	15 min
Adenovirus*	1.00%	5 min
Papovavirus/Polyomavirus/SV40*	0.50%	1 h
Poliovirus*	2.00%	15 min
Control of epidemics according to § 18 law for the protection against infections (IfSG)		
Disinfection of surfaces	4.00%	4 h
EURO-Norms (slight and strong contamination)		
EN 14476 (MNV)	1.00%	1 h
	2.00%	15 min
*DVV 2015 (Susp. Test)		

Important!

Place objects to be disinfected in the solution in such a way that all surfaces are in contact with the product. Close container tightly with lid. **Rinse instruments with water and endoscopes with germ-free water after disinfection /cleaning.**

Prepare dilution shortly before use (as recommended in the directives for Hospital Hygiene and Infection Prevention of the Robert-Koch-Institute, Berlin).

The solutions ready for use are efficient for at least **fourteen days** provided they are free of organic residue (disinfection of pre-cleaned instruments).

Shelf-life of the non-contaminated ready-for-use solutions when stored adequately (in a cool and dry place) and in clean containers with a tightly closed lid: 4 weeks.

Lysoform – Hygiene with System

concentration dosing system for 8 l of water	0.25 %	0.5%	1.0%
dosing pump (1 stroke \cong 20 ml)	1 stroke	2 strokes	4 strokes
Lysoform dosing devices	Adjust once to obtain the desired concentration		

Important: Always fill container with tap water first (at room temperature), then add disinfectant. The amount of water is calculated by subtracting the amount of disinfectant to be used from the total amount of solution wished for.

WARNINGS / STORAGE

- Keep disinfectants in a cool but frost free, well ventilated and dry place and out of reach of children.
- Wear suitable protective clothing, adequate protective gloves and eye protection when handling the product
- Do not use beyond expiry date.
- Shelf life of the solution ready for use free from organic substances: 14 days.
- Air the room after using as a surface disinfectant; Use in tightly closed containers (instrument disinfection).

PROPERTIES

- Bactericidal, tuberculocidal, mycobactericidal and virucidal,
- Broad spectrum of activity and good compatibility with surfaces, equipment, textiles, instruments
- Good cleaning power without addition of detergents, even in presence of organic substances like lipids, albumin, serum
- Freshly scented
- Resists hard water
- Biodegradable

SPECIAL PRECAUTIONS

a) Incompatibilities:

Adding detergents to **Aldasan® 2000/Aldosan® 2000-CH** can impair or neutralize its disinfecting properties.

b) Counter-measures in case of improper use: Should **Aldasan® 2000/Aldosan® 2000-CH** come into contact with the mucous membranes (eyes, throat, stomach etc.) seek medical advice immediately. Dilute product without delay.

If swallowed, drink **plenty** of water; after contact with the eyes, rinse immediately with **plenty** of water.

c) Miscellaneous:

After completion of the disinfection procedure rinse objects/surfaces which come into contact with food with drinking water.

d) Use **Aldasan® 2000/Aldosan® 2000-CH** safely. Always read the label and product information before use.

Certified and mentioned in the relevant lists.

STANDARD PACKAGING / TRANSPORT UNITS

6 x 1 litre bottles

2 x 5 litre cans

1 x 20 litre can

Dosing pumps for 5 l cans (1 stroke \cong 20 ml)

Dosing device DH 91

Article number

6241

6255

6260

3135

9174

Lysoform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Str. 133, D – 12247 Berlin
Tel.: +49-30-77992-0 Telefax: +49-30-77992-219
www.lysoform.de, e-mail: kontakt@lysoform.de

0418-16